



Enterprise Audit Trail

User Manual

Software Version 2.2

6006972-8EN
Revision 1
US English

General User Documentation

© 2026 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.
GE HealthCare proprietary. Reproduction and/or distribution is prohibited.

Enterprise Audit Trail

User Manual

6006972-8EN

Revision 1

February 2026

Proprietary to GE HealthCare - Class A

This material is the property of GE HealthCare and contains information that is restricted to use by GE HealthCare personnel. Any reproduction in whole or in part is strictly prohibited.

GE HealthCare Confidential. Disclosure to non-GE HealthCare parties is strictly prohibited.



Manufacturer:

GE Healthcare

500 W. Monroe Street

Chicago, IL 60661 USA

Revision History

Revision	Date	Reason for Change
1	February 2026	Initial Release.

Contents

1 Legal Information	5
1.1 Copyright Notice.....	5
1.2 Confidentiality	5
1.3 Trademark Notes	5
2 About this Manual.....	6
2.1 Intended Audience	6
2.2 User Interface Note	6
2.3 Demographic Information Disclaimer	6
2.4 Safety.....	6
2.4.1 Terminology	6
2.5 General Definition of Symbols.....	7
2.6 Access the User Manual Online.....	7
3 About Enterprise Audit Trail	9
3.1 Product Description.....	9
3.2 Minimum Screen Resolution	9
3.3 Browser Compatibility.....	9
3.4 Privacy and Security Considerations	10
3.5 Terms and Abbreviations.....	10
4 Getting Started.....	11
4.1 Signing in to the Application.....	11
4.2 Resetting the Login Password	13
4.3 Signing Out of the Application.....	14
5 Navigating the Home Page	15
5.1 Home Screen Menu Bar	15
5.2 Side Navigations	16
5.3 Footer	16
6 Using the Audit Event Page	17
6.1 Searching for Audit Events	18
6.1.1 Searching Events by Date Range.....	18
6.1.2 Searching Events by User and Patient Details	19
6.1.3 Searching Events by Activity Details	19
6.2 Viewing Audit Events.....	21
6.2.1 Hiding and Unhiding the Columns	22

6.2.2 Viewing User Details of Audit Event	22
6.2.2 Viewing Additional Details of Audit Event.....	23
6.2.3 Viewing Raw Audit Messages.....	23
6.3 Exporting Audit Logs.....	24
7 Troubleshooting.....	25

1 Legal Information

1.1 Copyright Notice

All Licensed Software and Documentation are protected by the copyright laws of the United States and by applicable international treaties.

© 2026 GE HealthCare Company. All rights reserved.

1.2 Confidentiality

This document is the confidential property of GE HealthCare and/or its affiliated entities. It is furnished to, and may only be used by, customers and their employees under a written agreement with GE HealthCare and may only be used in accordance with the terms of that agreement. The access and use of this document is restricted to customers and their employees. The user of this document agrees to protect the confidentiality of the information contained herein and GE HealthCare's proprietary rights as expressed herein and not permit access to this document by any person for any purpose other than as an aid in the use of the GE HealthCare software. In no case may this document or any portion hereof be accessed, made available, examined, or copied for the purpose of developing, marketing, or supporting any system or computer program similar to the GE HealthCare software. No part of this document may be copied without the prior written permission of GE HealthCare. The information in this document is subject to change by GE HealthCare without notice.

1.3 Trademark Notes

GE and the GE Monogram are trademarks of GE HealthCare Company. DoseWatch is a trademark of GE HealthCare Company or one of its subsidiaries. Windows, the Windows logo, Internet Explorer, and Microsoft Edge are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries. Android and Google Chrome are registered trademarks or trademarks of Google, LLC. macOS and Safari are trademarks of Apple Inc., registered in the U.S. and other countries. IOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license by Apple Inc. Firefox is a trademark of the Mozilla Foundation in the U.S. and other countries.

All other trademarks are the property of their respective owners.

2 About this Manual

2.1 Intended Audience

The Enterprise Audit Trail User Manual is intended for users and consumers of audit data, such as Legal or IT teams, Department or Service Line Managers, COOs, Process Improvement personnel, and, in some cases, Technologists and Medical Physicists.

2.2 User Interface Note

The screen captures in this document are intended for demonstration purposes only and may not always be fully identical to the actual user interface (e.g., in color scheme). Use these screen captures as guides.

2.3 Demographic Information Disclaimer

Any patient identifiable information contained in this manual is fictitious information generated for purposes of illustration only (including, but not limited to name, age, gender, date of birth, medical record number and exam date).

2.4 Safety

2.4.1 Terminology

The terms *Warning* and *Caution* are used throughout this manual to point out hazards and to designate a degree or level of seriousness. *Hazard* is defined as a source of potential injury to person.

Notice type	Description
Warning	Indicates a potentially hazardous situation, which, if not avoided, could result in death or serious injury.
Caution	Indicates a potentially hazardous situation, which, if not avoided, may result in minor or moderate injury.
Important	Indicates important information pertaining to a non-hazardous situation.
Note	Indicates informational messages providing application tips or other useful information that do not pertain to a hazard.

2.5 General Definition of Symbols

Symbol	Definition
	This symbol indicates the AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY of the product.
	This symbol indicates the MANUFACTURER of the product.
	This symbol indicates the DATE OF MANUFACTURE of the product.
	This symbol indicates that the operator should CONSULT INSTRUCTIONS FOR USE for further information.
	This symbol indicates the manufacturers BATCH CODE , or lot number of the product.
	The symbol for CAUTION highlights the fact that there are specific warnings or precautions associated with the device.
	The product bears CE mark indicating its conformity with the provisions of the Council Directive 93/42/EEC, concerning medical devices and fulfills the essential requirements of Annex I of this directive.
	The WEEE symbol indicates that the waste of electrical and electronic equipment must not be disposed of as unsorted municipal waste and must be collected separately.
Rx Only	For products distributed in the US, the symbol for “ Rx Only ” indicates: Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

2.6 Access the User Manual Online

This user manual is also available at GE HealthCare *Customer Documentation Portal* (*Manuals & Documents* section): <https://www.gehealthcare.com/support/manuals>.

You can order a paper copy of this user manual at no additional cost. Please, send a request for a copy to your sales or service representative. In application of the *EU Commission Regulation on electronic instructions for use of medical devices*, in the European Union, you can expect your request to be handled within 7 days.

1. Open *Customer Documentation Portal* (*Manuals & Documents*):
<https://www.gehealthcare.com/support/manuals>
2. Enter the part number of this user manual in the search field (**Enter product info** field).
 - 2.1. Optionally, select a product (or products) from the **Browse by product** field.
 - 2.2. Optionally, select a language (or languages) from the language selector (the default selected language on the drop-down field **English**).
3. Click the search icon.

 GE HealthCare

ProductsServicesResourcesSpecialtiesInsightsShopAbout Us



Manuals & Documents

 Support  Sign Out



 Browse by product

English 

All 13566

Service 3911

User/Operation 6306

Installation 2393 

Others 956 



NOTE

Alternatively, the user manual is also available for download within the Enterprise Audit Trail application. Click the  **Help** link on the side navigation of the Home page to open or download this user manual.

3 About Enterprise Audit Trail

3.1 Product Description

Within the GE Healthcare ecosystem, the Enterprise Audit Trail (EAT) service functions as a centralized single source of truth for application and clinical audit events, eliminating the need to maintain multiple audit repositories. The Audit Trail Repository supplies the messaging infrastructure that enables IHE-compliant actors and applications to securely receive, store, and retrieve audit messages. It provides features commonly required by IHE-compliant products such as message queues, a persistent data store, message search, and reporting.

The EAT service is scalable, fault-tolerant, and secure, with the ability to handle large volumes of audit events and provides insights into system activity. EAT receives and stores audit events as Fast Healthcare Interoperability Resources (FHIR) audit logs. For display, EAT reformats the events in a more user-friendly manner.

For more information on IHE-compliant products and FHIR audit logs, search for *IHE audit log guidelines*.

3.2 Minimum Screen Resolution

The minimum recommended resolution for client workstations to access Enterprise Audit Trail is 1280 x 1024. This resolution ensures that the data grid will not show horizontal or vertical scrollbars for the default page size. Horizontal or vertical scrollbars may be displayed for any other screen resolution or page size.

3.3 Browser Compatibility

The following browsers are used to run Enterprise Audit Trail:

- Computers with Windows® operating system
 - Microsoft Edge® and Internet Explorer® 11



NOTE

Compatibility mode for Internet Explorer® is not supported.

- Google Chrome™.
- Computers with macOS® operating system
 - Google Chrome™.
 - Firefox® and Firefox ESR
 - Safari®



NOTE

Google Chrome™ is the preferred browser.

3.4 Privacy and Security Considerations

The Enterprise Audit Trail receives ATNA-compliant (FHIR feed) audit log messages generated by other IHE products/components and stores them securely in a permanent data store. The content of the messages may contain personal information (PI) such as patient IDs (MRN, SSO), patient name, user IDs, and names. The Enterprise Audit Trail collects and stores this PI in its own data repository (Healthlake Datastore) for future use. The Enterprise Audit Trail (database) is hosted in AWS Cloud environments and managed within GEHC.

To ensure secure communication and to maintain data integrity of the audit log messages, the Enterprise Audit Trail uses a secured HTTPS protocol. The security risk assessment is evaluated based on the intended use of this component within a secure network boundary.

The searching and reporting of the audit messages are performed using the web-based user interface of the Enterprise Audit Trail tool. The user authentication to grant access to the EAT Web UI is managed by IDAM authentication services. A typical user of the EAT Web UI is an IT person or a supervisor who queries the log messages for auditing purposes.

The FHIR audit log repository enables the user to securely store the audit log events after validation, provides capabilities to do filtering, retrieval, querying and searching the events based on selected fields such as Date, Event action, etc.

3.5 Terms and Abbreviations

Term	Definition
EAT	Enterprise Audit Trail
ATNA	Audit Trail and Node Authentication (IHE profile)
FHIR	Fast Healthcare Interoperability Resources
IHE	Integrating the Healthcare Enterprise (initiative in the form of a standardization body that aims at promoting interoperability in healthcare systems)
HTTPS	Hypertext Transfer Protocol (Secure)
AWS	Amazon Web Service
IDAM	Identity and Access Manager
UI	User Interface
URL	Uniform Resource Locator
PI	Personal Information

4 Getting Started

4.1 Signing in to the Application

Consumers are required to log in with their credentials. The EAT UI is integrated with IDAM, and users with the appropriate subscription will be able to access the application upon authentication. Use this procedure to sign in and out of the EAT application.



NOTE

By default, no users have access to the EAT application. To request access or remove access, raise a support ticket with GEHC.

1. Type the URL in your web browser and press **Enter** to start the EAT application. The EAT sign-in page appears.

2. Type your username, an 'at' (@) symbol and the domain name of the organization with (with no spaces) in the **Username + '@' + organization domain** field.
That is, <username>@<domain name of the organization>.
3. Click **Continue**.

In the **Username + '@' + organization domain** field is auto filled with your username.



NOTE

The organization domain after the '@' symbol is not editable at this step, only the username can be edited. To change the organization domain, click the **Change**

organization link below the **Sign in** button.

4. In the **Password** field, type your password.
5. Click **Sign In**.

**NOTE**

If you have entered your username and/or password incorrectly or your user account has been deleted since your latest session for any reason, the application displays the following error message:

```
"Wrong username and/or password"
```

6. Retype your correct username and password.
7. Click **Sign In**.

If the problem persists, but you are confident that you are using the correct username and password, contact your administrator for assistance.

**NOTE**

You have five attempts to enter your password. If you enter your password incorrectly five times, your account is locked for 5 minutes after the first five failed sign-in attempts. You can retry signing in after your account has been unlocked. The account lock occurs after every five failed sign-in attempts, but the waiting time for the locked account to be unlocked is doubled after every account lock. The duration is five minutes after the first account lock, 10 minutes after the second one, and so on.

If you forgot your password, click the **Forgot password** link below the **Sign in** button to reset your password. For detailed instructions, see *GE HealthCare Identity and Access Manager User Manual*.

**NOTE**

By signing in, you accept the GE HealthCare IDAM **Terms of use**, **Privacy Policy**, and **CA Privacy Policy**. Click the **Terms of use** and/or the **Privacy Policy** and/or the **CA Privacy Policy** link to see more details.

**NOTE**

GE HealthCare IDAM uses cookies to improve browsing experience. Click the **Cookie policy** link to see more details.

**NOTE**

Click the **Get Help** link to open or download the GE HealthCare Identity and Access Manager user manual.

8. On the Authentication screen, complete your login authentication through the configured identity provider (IDP) or through e-mail based one-time password.


On your first sign-in to the application, the Consents for Identity and Access Manager application page opens, displaying the list of profile information the application requests or requires to work in the intended way. The **Allow all** check box is selected by default, so all optional information attributes are preselected.

9. Click **Submit** to sign in to the application.

4.2 Resetting the Login Password

Use the following procedure to reset the login password:

1. Type the URL in your web browser and press **Enter** to start the EAT application. The EAT sign-in page appears.
2. Type your username, an 'at' (@) symbol and the domain name of the organization with (with no spaces) in the **Username + '@' + organization domain** field and then click **Continue**. The IDAM Sign-in page appears.
3. Click the **Forgot password** link on the sign-in page. The Forgot password page appears.



4. The **Username + '@' + organization domain** field is auto-filled with your username if you have provided your details correctly before selecting the **Forgot password** link.
If your username and domain are not automatically filled in, type your username + '@' + organization domain.
5. In the **Captcha** field, type the characters exactly as they appear above the field (uppercase/lowercase, numbers).
If one or more characters appear illegible to you, click the **Refresh captcha** link.
6. Click **Submit**. An email is sent to the registered email ID.
7. Click the **Reset password** link in the email you have received. The **Reset password** page appears.
8. Type your new password in the **New password** field.
9. Type your new password again in the **Confirm password** field.

10. Click **Submit**.

After you reset your password successfully, sign-in using your new password.

4.3 Signing Out of the Application

To sign out of the application:


1. On the application **Home** page, click the name drop-down to display the options.
2. Click **Sign Out**.

The EAT sign-in page appears.

5 Navigating the Home Page

The EAT **Home** page is the landing page of the application after you successfully log in. Use this procedure to access the EAT **Home** page:

1. Sign in to the EAT application. For additional information, see [Signing in to the Application](#).
2. The EAT **Home** page opens.




The following features are available on the EAT **Home** page:

- [Home Screen Menu Bar](#)
- [Side Navigations](#)
- [Footer](#)


5.1 Home Screen Menu Bar

Following are the basic navigations and features available on the EAT Menu bar:

On the upper right corner of the Home page, click the  **<username>** profile drop-down option to perform the following actions:



- Change appearance between light and dark mode – Click  **Light** or  **Dark** to change the appearance of the application accordingly.
- View the **About** page – Click the  **About** option to view the version, build number, and manufacture details of the EAT application.



- **Sign out** – Click the  **Sign Out** option to sign-out from the application.

5.2 Side Navigations

Following are the basic navigations and features available on the side navigation:

- Switch between the Home Page and Audit Events page – Click either  **Home** or  **Audit Event** icon to open the respective windows.
- **Help** – Click the  **Help** icon to open the EAT user manual.

5.3 Footer

Following links are available at the footer of the EAT application:

- **Privacy** – Click the **Privacy** link to open the GE HealthCare's Privacy Policy page to understand how your personal data is collected, used, stored, and protected by the organization.
- **Terms of use** – Click the **Terms of use** link to open the GE HealthCare's Terms and Conditions page to understand the rules, responsibilities, and limitations associated with using the organization's applications or services.


6 Using the Audit Event Page

The Audit Event window provides an interface to query the audit log data stored in the EAT database and display the results.

The Audit Event window provides the following options:

- Search audit messages stored in the EAT database based on search criteria.
- Display the audit messages searched by the user.

To access the Audit Event window, click the  **Audit Event** icon in the side navigations. The Audit event window opens.



Audit event
Track and review important activities in one place

Find an activity

Date range
User and patient details
Activity details

Reset Find

Activities Found (58 records found)
3 column(s) hidden Unhide

Activity category	User ID	User name	Action taken	Patient name	Date & time	Details	Status
Login			Execute		2026-01-04T18:04:...		Success
Login			Execute		2026-01-04T15:01:...		Success
Login			Execute		2026-01-01T16:14:...		Success
Logout			Execute		2026-01-01T16:12:...		Success
Login			Execute		2026-01-01T16:11:...		Success
Login			Execute		2025-12-26T13:34:...		Success
Login			Execute		2025-12-26T13:31:...		Success
Login			Execute		2025-12-26T10:51:...		Success

The Audit Event window comprises two sections: **Find an activity** and **Activities Found**.

- **Find an activity** – This section helps a user to search for an audit event. For more information, see [Searching for Audit Events](#).
- **Activities Found** – This section displays the audit events as per the search criteria. If the search activity is not performed, this section displays all of the audit events recorded in the application. For more information, see [Viewing Audit Events](#).



NOTE

The audit events displayed in the **Activities Found** table are limited to the past six months.

6.1 Searching for Audit Events

EAT supports the following search capabilities:


- [Searching Events by Date Range](#)
- [Searching Events by User and Patient Details](#)
- [Searching Events by Activity Details](#)

6.1.1 Searching Events by Date Range

The most basic search for audit events is performed by selecting a date range.

To search events by a date range:

1. On the **Audit events** window, under the **Find an activity** section, expand the **Date range** section.



2. Searching is limited to the past six months. To search for all data up to the present, leave the **To** and **To time** fields blank. Enter the search criteria as follows:
 - **From** – Enter start date using MM/DD/YYYY format.
 - **From Time** – (optional) Enter start time using 24-hour clock format HH:MM:SS.
 - **To** – Enter end date using MM/DD/YYYY format.
 - **To Time** – (optional) Enter end time using 24-hour clock format HH:MM:SS.

NOTE




Searching is limited to the past six months. If you want to view the data beyond six months, raise a support request to GEHC.

3. Click **Find** to perform the search. The search result will be displayed in the **Activities Found** table.
4. Click **Reset** to clear the search criteria values and results.

6.1.2 Searching Events by User and Patient Details

Follow the process below to search for the events by User and Patient Details:

1. On the **Audit events** window, under the **Find an activity** section, expand the **User and patient details** section.




2. Enter the details in the fields as per the following:
 - **User ID** – Enter the ID of the person who was involved in an action or event. It points to the person, device, or system that performed or participated in the activity being audited.
 - **Patient ID** – Enter the ID through which patient was involved in an action or event.
 - **User name** – Enter the name of the person who was involved in an action or event. It points to the person, device, or system that performed or participated in the activity being audited.
 - **Patient name** – Enter the name through which patient was involved in an action or event.
3. Click **Find** to perform the search. The search result will be displayed in the **Activities Found** table.
4. Click **Reset** to clear the search criteria values and results.

6.1.3 Searching Events by Activity Details

Follow the process below to search for the events by activity details:

1. On the **Audit events** window, under the **Find an activity** section, click the drop-down

button next to **Activity details** section to expand the section.



The screenshot shows a search interface titled "Find an activity" with a clock icon and an information icon. It features three expandable sections: "Date range" (collapsed), "User and patient details" (collapsed), and "Activity details" (expanded). Below these sections are two input fields: "Source of activity" with a placeholder "Enter Source of activity" and "Item ID" with a placeholder "Enter Item ID". A note below the "Item ID" field states "Support Study ID, Series ID and Accession number". At the bottom are "Reset" and "Find" buttons.

2. Enter the details in the fields as per the following:
 - **Source of activity** – Enter the application or system name that sends the audit log of audit events to search (For examples of valid activity sources, note the App ID included in the audit event details). This allows for an App ID based search.
 - **Item ID** – Enter alphanumeric text corresponding to the specific event such as a study ID, series ID, accession number, patient record, lab result, or document number.
3. Click **Find** to perform the search. The search result will be displayed in the **Activities Found** table.
4. Click **Reset** to clear the search criteria values and results.

6.2 Viewing Audit Events

The **Activities Found** table displays the list of audit events or the audit events that comply with the search criteria. Audit events are displayed in reverse chronological order (newer events in the top rows and older events in the bottom rows.)

The total number of audit events is displayed next to **Activities Found**. By default, only partial information is displayed in each search results row. To display additional information, hover the mouse over a row, then click the arrow next to each row.

Activities Found (59 records found)										
1	2	3	4	5	6	7	8	9	10	11
ID	Event Time	Activity category	User ID	User name	Action taken	Patient ID	Patient name	Date & time	Details	Status
▼		Login			Execute			2026-01-04T22:14:...		Success
▼		Login			Execute			2026-01-04T18:04:...		Success
▼		Login			Execute			2026-01-04T15:01:...		Success
▼		Login			Execute			2026-01-01T16:14:...		Success
▼		Logout			Execute			2026-01-01T16:12:...		Success
▼		Login			Execute			2026-01-01T16:11:...		Success
▼		Login			Execute			2025-12-26T13:34:...		Success
^		Login			Execute			2025-12-26T13:31:...		Success

The following table describes the different columns in the **Activity Found** table.


#	Columns	Description
1	ID	ID of the event
2	Event Time	The local time of the location when the actual event occurred
3	Activity category	Category of the activity
4	User ID	ID of the person who performed the action
5	User name	Name of the person who performed the action
6	Action taken	Indicator for type of action performed during the event that generated the event, for example: - Create - Create a new database object, such as placing an order - Read/View/Print - Display or print data, such as a doctor census - Update - Update data, such as revise patient information - Delete - Delete items, such as a doctor master file record - Execute - Perform a system or application function such as log-on, program execution or use of an object's method, or perform a query/search operation.
7	Patient ID	ID of the patient
8	Patient name	Name of the patient
9	Date & time	Date and time when the audit record was created
10	Details	Details of the raw audit message
11	Status	Status of the event

6.2.1 Hiding and Unhiding the Columns

You can use filters to hide or unhide the columns in the Activities Found table as per your preference.

To hide or unhide the columns:

1. On the top right of the Activities Found table, click the  icon. The **Hide Columns** pop-up opens.



2. Select the column heading checkboxes that you want to hide from the **Activities Found** table or unselect the headings that you want to display.
3. Click **Save**. The **Activities Found** table is displayed as per your settings.



NOTE

If you have chosen to hide the columns, the **Unhide** link appears on top of the table. Click **Unhide** to unhide all the columns. Click the  icon to revert your changes to your previous selection.

6.2.2 Viewing User Details of Audit Event

It is sometimes necessary to view details of the users involved in an Audit Event. Some activities involved multiple users.

To view User details:

1. In the **Activities Found** table, click the User ID hyperlink of an audit. The **Users** pop-up opens and displays the number of user IDs involved in the audit event.

Users

User ID

User Name

Type

Roles

Network

human user

AuditAdmin

unknown

Close

2. Click **Close** to close the pop-up window.

6.2.2 Viewing Additional Details of Audit Event

You can see the additional details of each audit event either by clicking the active links or by expanding each audit event. This enables you to view details of multiple values of users, IDs, or systems linked to the audit event. For example, there could be multiple User IDs since multiple actors can be involved in an audit event.

To view the additional details of an audit event:

1. To view additional details of an activity, expand the activity. The individual audit events are displayed.

Activity categ...	User ID	User name	Action taken	Patient name	Date & time	Details	Status
^ Query			Execute		2026-02-16T0...		Success
Item ID	Activity type		Source of activity		Purpose of use		
-	search		Device/eatclientid		-		
^ Query			Execute		2026-02-16T0...		Success
Item ID	Activity type		Source of activity		Purpose of use		
-	search		Device/eatclientid		-		


2. Click **Close** to close the pop-up.

6.2.3 Viewing Raw Audit Messages

The EAT application enables you to view the raw FHIR R4 audit message.

To view the raw audit message:

1. Click the  icon under the **Details** column of an audit event. The **Details** pop-up opens.



2. Close the pop-up to return to the Audit events page.

6.3 Exporting Audit Logs

To obtain an electronic version of the audit logs, raise a support request to GEHC. This version of EAT doesn't support exporting or saving audit logs due to patient privacy restrictions.

7 Troubleshooting

1. Invalid username or password

If the username or password entered on the *IDAM Sign-in* page is invalid or not found, the system displays the following message:

Error Message:

Login failed!

Please check the username and password and try again.

Possible resolutions:

- If the username is correct, you can reset the password at the password entry step on the Sign-in screen. For instructions, see [Resetting the Login Password](#).
- If you do not know the correct username to use, contact your **Organization Administrator** for assistance who can verify the account or force a password reset.

2. Unrecognized username domain

During any login attempt to or through IDAM, the system first displays a page for selecting your organization. This is the page that asks for the login information in the following format: `Username + @ + organization domain`. If you have entered the customer domain incorrectly, then the system cannot find the organization domain and displays the following message:

Error Message:

Authentication error!

Something went wrong during the authentication process. Please try signing in again.

Possible resolutions:

- If the organization is configured correctly, it is possible that you are trying to use the wrong user account. If you do not know the correct username to use, contact your **Organization Administrator** for assistance.
- In some cases, the browser's existing session cookies can also cause such authentication errors. You can resolve this issue if you clear your browser's cache and cookies, or close all open browser windows, then open a private browser window, and try again.

3. Connection problem

The following error appears when there is a network issue:

Error Message:

Please check your network and try again.

Possible resolutions:

- Check your internet connection.
- The server may be temporarily unavailable.

- Your firewall or proxy may be blocking the request.

4. **Invalid input**

Error Message:

Please check your details and try again.

Possible resolutions:

- Tenant ID is missing or invalid, contact your Administrator.

5. **Login failed**

Error Message:

Please check your credentials and try again.

Possible resolutions:

- Your session may have expired - try refreshing the page.
- Invalid API credentials. You may need to log in again.

6. **Access denied**

Error Message:

You don't have permission to access this resource.

Possible resolutions:

- You don't have permission to access this data.
- Your account may not have the required privileges.

7. **No data found**

Error Message:

It may have been moved or deleted.

Possible resolutions:

- The requested audit events may not exist.
- Check if the service is properly configured.

8. **Too many requests**

Error Message:

Please wait a moment and try again.

Possible resolutions:

- You've exceeded the rate limit - wait for some time before you try again.
- Try reducing the frequency of your searches.
- The system is temporarily limiting requests.

9. **Server Error**

Error Message:

Please try again later.

Possible resolutions:

- The search query might have failed.
- There's a problem with the server.

10. Service unavailable right now

Error Message:

Please try again in a few minutes.

Possible resolutions:

- The server is temporarily down for maintenance.
- Network connectivity issues between services.

11. Service overloaded

Error Message:

Please try again later.

Possible resolutions:

- The service is temporarily overloaded.
- Server maintenance is in progress.

12. Request timed out

Error Message:

Please check your network and try again.

Possible resolutions:

- The server took too long to respond.
- Try simplifying your search criteria.
- Network connectivity may be slow.



www.gehealthcare.com